

# Long-term clinical evidence with Ankylos<sup>®</sup> implants

NEWLY PUBLISHED RESULTS WITH DENTSPLY SIRONA IMPLANTS PRODUCTS

Up to **17** years  
OF CLINICAL RELIABILITY

**18 945**  
ANKYLOS IMPLANTS

A retrospective analysis of Ankylos implants placed in a large patient population showed successful long-term clinical results when using a flapless surgical approach.



**7783**  
patients

Implant survival  
rate after  
**1 year** - 98.5%  
**5 years** - 96.7%  
**10 years** - 93.0%



Minimally  
invasive  
surgery



Jesch P, Jesch W, Bruckmoser E, et al. An up to 17-year follow-up retrospective analysis of a minimally invasive, flapless approach: 18 945 implants in 7783 patients. Clin Implant Dent Relat Res 2018;E-pub Feb 15, doi: 10.1111/cid.12593. [Abstract](#)



[www.dentsplysirona.com/implants/science](http://www.dentsplysirona.com/implants/science)

**Dentsply  
Sirona**  
Implants